

HIFM+RF Beauty muscle instrument

Muscles make up about 35% of the body, yet most slimming devices on the market focus only on fat, not muscle. Traditionally, only injections and surgery have been used to improve the shape of the buttocks. In contrast, the HIFM Beauty Muscle Instrument utilizes advanced High-Intensity Focused Magnetic (HIFM) technology to stimulate motor neurons, causing muscles to expand and contract in a way that cannot be achieved through regular exercise. A single 30-minute session delivers energy pulses that induce 30,000 intense muscle contractions, helping fat cells metabolize and break down effectively. Alongside muscle strengthening, it offers a cutting-edge experience in body shaping. The device has received FDA and CE certifications, and its safety and effectiveness are widely recognized.

FEATURES

1. New high-intensity focused magnetic vibration + focused RF; improves postpartum rectus abdominis separation.
2. Adjustable settings offer different muscle training modes.
3. Dual handles support synchronized operation for two users simultaneously, suitable for both men and women.
4. HIFM and RF energy can be applied independently or combined, allowing flexible treatment with one or both energy types; RF dual-channel control enables independent adjustment of treatment parameters.
5. RF heat energy is released from the inside out, ensuring no damage to skin or muscles, with a warm and comfortable treatment experience.
6. Safe and non-invasive, with no current, hyperthermia, or radiation, and no recovery period required.
7. Time-efficient: 30 minutes of treatment induces 30,000 muscle contractions (equivalent to 30,000 sit-ups or squats).
8. Extensive studies confirm remarkable results; only 4 treatments within two weeks are needed to see noticeable contouring effects.
9. An air-cooling system prevents overheating of the treatment head, enabling continuous operation, enhancing durability, and ensuring stable energy output.



Specifications

Product name	HIFM+RF Beauty muscle instrument		
Magneti vibration intensity	7 Tesla		
RF temperature	40-50°C	RF frequency	40.68M
Input voltage	AC110V/230V		
Output power	300W-3000W		
Output frequency of Magnetic resonance wave	3-150HZ		
Fuse	20A		
Host size/weight	52×39×34cm		
Size of flight shipping case	47×65×63cm		
Gross weight	39kg		

6D LIPO LASER

Laser Lipo emits low levels of laser energy, which creates a chemical signal in the fat cells, breaking down the stored triglycerides into free fatty acids and glycerol and releasing them through channels in the cell membranes. The fatty acids and glycerol are then transported around the body to the tissues that will use them during metabolism to create energy. This process of fatty acid release is a natural response of the body when the body needs to use stored energy reserves; thus, Lipo is not creating any unnatural reaction in the body nor does it affect or damage any surrounding structures such as skin, blood vessels and peripheral nerves. A period of exercise post-treatment will ensure the complete metabolism and elimination of the freed fatty acids from the body.

Features

1. The 6-D Lipo Laser is equipped with 230 diodes from Japan Mitsubishi , each diode with a output of 209mw.
2. A laser diode that is able to generate light with a wavelength of 635nm, 650nm, 780nm, 810nm, 940nm, 980nm . 10 laser pads (8 large and 2small pads) designed to work on different parts of the body.
3. The device is fully computerized and equipped with an intelligent easy-to-use 14.1-inches touchscreen . Heavy duty laser unit design
4. 6-D Lipo laser offers additional rental options for sales agents around the world. IC card and remote control rental.



Specifications

Model	6-D Lipo Laser
Wavelength	635nm+650nm+780nm+810nm+940nm+980nm
Output Energy	209mw/diode laser
Number of big paddles	8
Nuber of small paddles	2
Screen	14.1 inches
Voltage	110V-220V
Cooling	Air cooling
Mode of operation	CW pulse/Pulse
Warranty	2 years